

apteryx  
IMAGING

---

xv web<sup>®</sup>  
CAPTURE

# System Requirements

## XVWEB CAPTURE

---

XVWeb Capture and its related tools, components, and plug-ins are designed for 32- and 64-bit Windows® 10, Windows Server® 2008 / 2008 R2, Windows Server® 2012, and Windows Server® 2016 operating systems including Remote Desktop Services (formerly Terminal Services) and Citrix®.

Due to the enterprise nature of XVWeb Capture, server requirements are directly related to the size and expected growth of the database and the expected number of users that will be simultaneously using the application. Please contact Apteryx Technical Support for specific system requirements for your site at 877-APTERYX (877-278-3799).

### Ideal System Recommendations

- 2.0 GHz or higher processor (faster is always better, especially for servers)
- 4GB of RAM (again, more is always better)
- 500 GB Hard Drive
- CD-R/DVDR drive
- 100/1000Mbps or higher Local Area Network (TCP/IP preferred)

### Minimum System Requirements

The system must conform to the minimum requirements specified by Microsoft® for the operating system installed.

## REVISION HISTORY

---

Revision Level	Revision Date	Description of Revision	Effective Date	Revision Author Initials	Approver Initials	Doc Change Task #
1.0	04/05/2018	Initial release	04/16/2018	MG	PW	788
1.1	01/27/2020	Removed Windows 7 from "designed for" statement.	01/27/2020	MG	PW	884
1.2	01/27/2020	Removed Windows Vista from "designed for" statement.	01/27/2020	MG	PW	884
1.3	05/06/2021	Updated general document design, brand colors and logo/cover graphics, updated company URL.	05/24/2021	DC	PW	906
1.4	09/15/2023	Removed Windows 8 from "designed for" statement.	09/15/2023	NK	AS	1915