

apteryx

IMAGING

DataGrabber



System Requirements

DATAGRABBER

DataGrabber, Apteryx's patented bridge to practice management software, is designed for 32- and 64-bit Windows® 8, Windows® 10, Windows Server® 2008 / 2008 R2, Windows Server® 2012, and Windows Server® 2016 operating systems including Remote Desktop Services (formerly Terminal Services) and Citrix®.

Depending on the PM application, an ODBC Data Source may need to be specified on each workstation. Please contact Apteryx Technical Support for specific support information at 877-APTERYX (877-278-3799).

Application Requirements

- Installed practice management (PM) solution
- Appropriate Apteryx DataDrill to extract patient information from PM database
- Third-party dental or medical imaging software

Minimum System Requirements

The system must conform to the minimum requirements specified by Microsoft® for the operating system installed.

REVISION HISTORY

Revision Level	Revision Date	Description of Revision	Effective Date	Revision Author Initials	Approver Initials	Doc Change Task #
1.0	08/01/2011	Initial release	08/01/2011	SN		
1.1	04/05/2018	Removed Windows XP and Windows Server 2003 from "is designed for" statement. Added Windows 8, Windows 10, Windows Server 2012 and Windows Server 2016 to statement.	04/16/2018	MG	PW	788
1.2	01/27/2020	Removed Windows Vista and Windows 7 from "is designed for" statement.	01/27/2020	MG	PW	884
1.3	05/06/2021	Updated general document design, brand colors and logo/cover graphics, updated company URL.	05/24/2021	DC	PW	884